

# Covid-19 project - Digital lab journal

## Experiments

Study number	Name/Description	Date
SN0001	Trondheim #1 w/7 primer sets	02.09.2020
SN0002	Control #2 w/7 primer sets	02.09.2020
SN0003	Homemade vs kit buffer	03.09.2020
SN0004	Old vs new homemade buffer	03.09.2020
SN0005	One-step PCR, test of Taq Native vs Platinum	11.09.2020
SN0006	One-step PCR w/1X and 2X template conc.	14.09.2020
SN0007	One-step PCR w/Sarbeco and diff. Taq times	17.09.2020
SN0008	One-step RT-PCR of all primers	17.09.2020
SN0009	Finding optimal temp. for all primer sets	21.09.2020
SN0010	Reducing extension time to 30 seconds	21.09.2020
SN0011	Weekly run: long-term storage of mastermix	Week 39-44
SN0012	Standard curves of diluted template	24.-30.09.2020
SN0013	Reducing RT and anneal/extension time	01.10.2020
SN0014	Comparing setup10 with setup14	02.10.2020
SN0015	Standard curve of new Trondheim RNA	07.10.2020
SN0016	Test of kit before sending to Rikshospitalet	09.10.2020
SN0017	Enzyme/RNA impacts post-storage in RT	16.10.2020
SN0018	Fresh vs stored mastermix	16.10.2020
SN0019	2 hour storage in either fridge or freezer	21.10.2020
SN0020	Nuclear (nuTH01) and mitochondrial (Mito) primer sets for IC	21.10.2020
SN0021	Test of nuTH01, Mito and 3 PSA primer sets	22.10.2020
SN0022	Reproducing test performed by Rikshospitalet (lower template conc.)	23.10.2020
SN0023	Radium vs Riks probe concentrations w/random hexamer	23.10.2020
SN0024	Test of nuTH01 and Mito in different buffers	28.10.2020
SN0025	Contamination check of sterile water + diff. DNA ratios in IC primers	29.10.2020
SN0026	Cont. check of sterile water + IC with temp. gradient and 2X SW10 DNA	30.10.2020
SN0027	Trondheim and Oslo RNA - different concentrations and kit mix	02.11.2020
SN0028	New conc. of #21 and MS2, 2X template, gradient and 20 vs 25 µL	03.11.2020
SN0029.1	New setup - 2X Oslo or Trondheim RNA w/#21 and/or MS2, 20 vs 25 µL	03.11.2020
SN0029.2	SN0029.1 repeat with 5X RNA and 25 µL samples	03.11.2020
SN0030	Weekly run: setup30, all Covid-19 primer sets and MS2 IC	Week 46-
SN0031	Setup30 with 1X Trondheim RNA and 3X Spytt RNA	11.11.2020
SN0032	1:5 MS2 dilution series	12.11.2020
SN0033	Patients positive for Sars-Cov2 in Oslo RNA plate	12.11.2020
SN0034	Testing Bhadra et al (2020) protocol for RT with Taq only	19.11.2020
SN0035	Finding the optimal Mg <sup>+</sup> concentration with our mix	25.11.2020
SN0036	Finding the opt. Mg <sup>+</sup> concentration with our one-step and two-step mix	
SN0037	Testing the new and old batch of MS2 and TM primer/probe sets	02.12.2020
SN0038	Testing different concentration of Mg <sup>+</sup> with KrasX2 primer/probe set	02.12.2020
SN0039	Trying MS2 internal control in Oslo RNA and Trondheim RNA	02.12.2020
SN0040		
SN0041		
SN0042		



# Covid-19 project - Digital lab journal

## Trondheim #1 w/7 primer sets

SN0001

02.09.2020

### Aim:


Run 1st RT-PCR with Trondheim #1 sample (cDNA) and 7 primer sets. All parameters are based on preliminary data from Per O. Ekstrøm.

### Mix:


Fragment: Sars CoV	PCR		Grad 50-63C	Date	02.09.2020	
PCR volum, µl	20		# of reactions		9	
	Working					
	solutions		Total volume 180 µl		Desired	
	concentration		Volum		concentration	
H <sub>2</sub> O			155,16			
10X Thermopol uten MgSO <sub>4</sub>	0 mM MgCl		18,00			
MgSO <sub>4</sub>	200 mM		1,80		2 mM	
Primer forward	100 µM		0,72	Used 2,2 µL of each	0,4 µM	
Primer reverse	100 µM		0,72		0,4 µM	
Probe	100 µM		0,72		0,4 µM	
dNTP	100 mM		0,72		400 µM	
cDNA	10 ng				0,5 ng/µl	
BSA	100 %		1,80		1 %	
Dave Taq native	37,5 U/µl		0,36		0,075 U/µl	
Mashup	7,5 U/µl				0,075 U/µl	
Fragment: Sars CoV 2	PCR			Date	02.09.2020	
PCR volum, µl	20		# of reactions		70	
	Working					
	solutions		Total volume 1400 µl		Desired	
	concentration		Volum		concentration	
H <sub>2</sub> O			1083,60			
10X Thermopol #1	0 mM MgCl		140,00			
MgSO <sub>4</sub>	200 mM		14,00		2 mM	
Primer forward	100 µM				0,4 µM	
Primer reverse	100 µM				0,4 µM	
Probe	100 µM				0,4 µM	
dNTP	100 mM		5,60		400 µM	
cDNA	7 ng		140,00		0,7 ng/µl	
BSA	100 %		14,00		1 %	
Taq Hot	37,5 U/µl		2,80		0,075 U/µl	
Mashup	7,5 U/µl		0,00		0 U/µl	

### Results:


A



B



**C**



**First RT-PCR. (A)** A lot of background noise can be seen in the image. Reasons for this could be bubbles (remember to spin down plate at 1000 rpm before RT-PCR) or pipetting errors. **(B)** 2nd run with less noise. **(C)** Box plot of control sample results from the RT-PCR.

### Conclusion:

Primer set #7 and #20 shows the most promising results, with low cycle values (Cq/Ct) and high fluorescence. Primer set #20 also seem more specific than most of the other sets.



## Control #2 w/7 primer sets

SN0002

02.09.2020

### Aim:

Check if PCR results improves with 10  $\mu$ L mix compared to 20  $\mu$ L.  
Used pre-made cDNA.


### Mix:

Same as in SN0001, with some adjustments in regards to number of samples.

### PCR protocol:

- 50-61°C temperature gradient
- 50 cycles

### Results:



### Conclusion:

Some values were not included in the box plot as they represented outliers (a lot of noise).



## Homemade vs kit buffer

SN0003

03.09.2020

### Aim:

See if homemade buffer is as good/better than a kit buffer.  
Used pre-made cDNA from Control #1 RNA and Trondheim RNA.

### Mix:

PCR run with Per, 03.09.2020. Tested 2 different volumes (1-2  $\mu$ L) of cDNA (after synthesis, template volume pipetted directly into wells) and homemade buffer vs. kit buffer. Samples were tested with primer #7 and #20, as previous experiments have shown that these primers gives best effect (steepest curves, little background noise).


Fragment: Sars CoV 2 #7/20	PCR			Date	03.09.2020	
PCR volum, $\mu$ l	20		# of reactions		50	
	Working					
	solutions		Total volume 1000 $\mu$ l		Desired	
	concentration		Volum		concentration	
H2O			862.00			
10X Thermopol #1		0 mM MgCl	100.00			
MgSO4	200	mM	10.00		2	mM
Primer forward	100	$\mu$ M	4.00		0.4	$\mu$ M
Primer reverse	100	$\mu$ M	4.00		0.4	$\mu$ M
Probe	100	$\mu$ M	4.00		0.4	$\mu$ M
dNTP	100	mM	4.00		400	$\mu$ M
cDNA	10	ng	0.00	Adjust to 1X or 2X	0	ng/ $\mu$ l
BSA	100	%	10.00		1	%
Dave Taq native	37.5	U/ $\mu$ l	2.00		0.075	U/ $\mu$ l
Mashup	7.5	U/ $\mu$ l	0.00		0	U/ $\mu$ l

### PCR protocol:


- Same PCR settings as in SN0002

### Results:


A



**B**



**C**



**Homemade vs kit buffer. (A)** Regression-fitted amplification curve from experiment. **(B)** Box plot of Cq values when different amounts of cDNA from Control #1 RNA was used with primer set #7 or #20, and with different types of buffers. **(C)** Same conditions as in (B), except for the template. Here, cDNA from "Trondheim RNA" was used (purified RNA, received from Dr. Magnar Bjørås).

## Conclusion:

- Homemade buffer gave better results in Control cDNA samples than kit buffer, but not in Trondheim samples.
- Overall, 2X template gave more specific values with lower Cq - not that surprising as more viral RNA gets detected earlier/better.
- We conclude that we can safely use the homemade 10X Thermopol buffer (with Mg<sup>+</sup>) for later experiments.
- I added MgSO<sub>4</sub> by mistake in all buffers, but it didn't seem to induce any notable differences.



## Old vs new homemade buffer

SN0004

03.09.2020

### Aim:

Per made a really good batch of homemade buffer, but it is almost empty - we want to check if the fresh one is just as good for usage (the pH meter was broken when the latter was made).

### Mix:

Used primer sets #7 or #20, in addition to 2 different template concentrations (1X and 2X).


Fragment: Sars CoV 2 #7/20	PCR			Date	03.09.2020	
PCR volum, $\mu$ l	20		# of reactions		25	
	Working					
	solutions		Total volume 500 $\mu$ l		Desired	
	concentration		Volum		concentration	
H2O			431.00 or 436			
10X Thermopol old or new	0 mM MgCl		50.00	Only 1 buffer w/Mg <sup>+</sup>		
MgSO4	200 mM		5.00 or 0		2 mM	
Primer forward	100 $\mu$ M		2.00		0.4 $\mu$ M	
Primer reverse	100 $\mu$ M		2.00		0.4 $\mu$ M	
Probe	100 $\mu$ M		2.00		0.4 $\mu$ M	
dNTP	100 mM		2.00		400 $\mu$ M	
cDNA	10 ng		0.00	Adjust to 1X or 2X	0 ng/ $\mu$ l	
BSA	100 %		5.00		1 %	
Dave Taq native	37.5 U/ $\mu$ l		1.00		0.075 U/ $\mu$ l	
Mashup	7.5 U/ $\mu$ l		0.00		0 U/ $\mu$ l	

### PCR protocol:


- Same PCR settings as in SN0002

### Results:

A



**B**



**Old vs new homemade buffer. (A)** Regression-fitted amplification curve from experiment. **(B)** Box plot from the experiment.

## Conclusion:

- The results were quite varying but it seems that the fresh buffer gives just as good results as the old (good) batch. We can continue with the new one when the old runs out.
- 2X template gives a bit better Cq, consistent with our previous findings (see SN0003).





## One-step PCR, test of Taq Native vs Platinum

SN0005

11.09.2020

### Aim:

To test the possibility of doing the whole reaction in one step.

### Mix:

Same as in SN0001, with some adjustments:


- RNA (w/MashUp RT) vs cDNA (w/o RT)  
1-step 2-step
- Native Taq (Nat) vs Platinum Taq (Plat)
- Primer set #7 vs #20
- 8 conditions in total, 12 samples per condition

### PCR protocol:


- Same PCR settings as in SN0002
- 10 min cDNA synthesis at the start of the PCR protocol

### Results:

A



B



**One-step PCR test with different Taq polymerases. (A)** Regression-fitted amplification curve from the experiment. **(B)** Box plot from the experiment.



### Conclusion:

- No large difference between native and platinum Taq in #20 samples.
- Native Taq shows up a bit later in #7 samples but in parallel with Platinum Taq.  
This should not be affected by the temperature gradient as the differences are too small.
- From these results, it looks like one-step RT-PCR with our components work!  
Should be repeated for verification.



## One-step PCR w/1X and 2X template conc.

SN0006

14.09.2020

### Aim:

Repeat SN0005 to verify that one-step RT-PCR works, and check whether 2X template concentration improves the results.

### Mix:

Same as in SN0005, with some adjustments:


- Either 1X or 2X template concentration

### PCR protocol:

- Same as in SN0005

### Results:

A



### Conclusion:

- Could try to reduce cDNA synthesis time from 10 min to 5 min.
- Did not have enough Plat Taq/Hot start - only RNA samples got this polymerase, all cDNA samples got Nat Taq (some Plat Taq got diluted (2  $\mu$ L + 1.9  $\mu$ L H<sub>2</sub>O)).
- In general, 2X works better than 1X template
- The box plots do not show any significant difference between Nat and Plat Taq
- RNA samples have lower Cq/Ct than cDNA samples (more copies)
- Primer #20 shows most specificity/least variance
- Row 1 had no template (negative control). In well H01 (grey box in box plot) some fluorescence was detected at ~35 Cq, which have been excluded in the figure.



## One-step PCR w/Sarbeco and diff. Taq times

SN0007

17.09.2020

### Aim:

As one-step RT-PCR worked so far, we want to test the settings/mix with Sarbeco, the primer set used in standard testings in Norway. In this experiment, we also try to reduce the cDNA synthesis time to 5 min and include both 5 and 10 min incubation in the experiment (sealing half of the plate).

### Mix:

Fragment: Sars CoV Sarb	PCR		Grad 50-63C	Date	17.09.2020
PCR volum, µl	20		# of reactions		50
	Working				
	solutions		Total volume 1000 µl		Desired
	concentration		Volum		concentration
H <sub>2</sub> O			802.00		
10X Thermopol uten MgSO <sub>4</sub>	0	mM MgCl	100.00		
MgSO <sub>4</sub>	200	mM	10.00		2
Primer forward	100	µM	4.00		0.4
Primer reverse	100	µM	4.00		0.4
Probe	100	µM	4.00		0.4
dNTP	100	mM	4.00		400
RNA/cDNA	10	ng	50.00		0.5
BSA	100	%	10.00		1
Dave Taq native	37.5	U/µl	2.00		0.075
Mashup	7.5	U/µl	10.00		0.075


**S.O.E.:** Samples with 10 min RT incubation time got twice as much Taq! Repeated experiment.

### PCR protocol:


- Same as in SN0005, except that 10 min samples were incubated at 42°C for 5 min prior to adding the 5 min samples (no enzyme loaded in 5 min samples at the start). Plate was partly sealed with tape (Clas Ohlson).

### Results:

A



**B**



**One-step RT-PCR with Sarbeco primer set and different Taq polymerase times. (A)** Box plot based on regression-fitted values. **(B)** Box plot from based on single threshold Cq values.

### Conclusion:

- Although the 10 min incubation samples have lower Cq values, the 5 min samples are more specific.
- We will use 5 min RT incubation time from now on.



## One-step RT-PCR of all primers

SN0008

17.09.2020

### Aim:

Test all primers with the latest set-up (including 5 min RT incubation)  
Used pre-mixed primer/probe mix from strip.

### Mix:

I loaded the primer/probe mix into wells before adding the mastermix.

### One-step RT-PCR


Sample			Date	17.09.2021
PCR volum, $\mu\text{L}$	20		# of reactions	100
Reagents	Working sol. conc.	Total volume 2000 $\mu\text{L}$ Volume		Desired conc.
H <sub>2</sub> O		1627.98		
10X Thermopol uten MgSO <sub>4</sub>	0	200.00		
MgSO <sub>4</sub>	200 mM	20.00		2 mM
Primer forward	100 $\mu\text{M}$	0.01	Used pre-mixed primer/probe mix	0.0004 $\mu\text{M}$
Primer reverse	100 $\mu\text{M}$	0.01		0.0004 $\mu\text{M}$
SYBR green	100 $\mu\text{M}$	0.01		0.0004 $\mu\text{M}$
dNTP	100 mM	8.00		400 $\mu\text{M}$
RNA	10 ng	100.00		0.5 ng/ $\mu\text{l}$
BSA	100 %	20.00		1 %
Dave Taq native	37.5 U/ $\mu\text{l}$	4.00		0.075 U/ $\mu\text{l}$
Mashup	7.5 U/ $\mu\text{l}$	20.00		0.075 U/ $\mu\text{l}$

### PCR protocol:


- Same as in SN0006, with only 5 min RT incubation time

### Results:

A



**B**



**One-step RT-PCR of all primers. (A)** Box plot based on regression-fitted values. **(B)** The graph shows optimal temperature for annealing during RT-PCR for each primer/probe set.

## Conclusion:

- We will continue to use the temperature that gives the best (i.e. steepest) curve for each primer (see continuation in SN0009).



## Finding optimal temp. for all primer sets

SN0009

21.09.2020

### Aim:

To compare all primer/probe sets with their respective optimal conditions.  
Optimal temperature are based on data from project SN0008.

### Mix:

I used the same mix as in SN0008 but with some adjustments:

#### One-step RT-PCR

Fragment:	Trondheim RNA	Date	21.09.2020
PCR volum, µL	20	Temp. Gradient 55-62°C	# of reactions 12

Reagents	Working sol. conc.		Total volume 240 µL		Desired conc.	
			Volume			
H <sub>2</sub> O			192.48			
10X Thermopol uten MgSO <sub>4</sub>	0		24.00			
MgSO <sub>4</sub>	200	mM	2.40		2	mM
Primer forward	100	µM	0.96	Or 2.88 µL from mix	0.4	µM
Primer reverse	100	µM	0.96		0.4	µM
SYBR green	100	µM	0.96		0.4	µM
dNTP	100	mM	0.96		400	µM
RNA/cDNA	10	ng	12.00	Add after blank	0.5	ng/µl
BSA	100	%	2.40		1	%
Dave Taq native	37.5	U/µl	0.48		0.075	U/µl
Mashup	7.5	U/µl	2.40		0.075	U/µl

### PCR protocol:

- Same as in SN0006, but with a temperature gradient fit for the different primer/probe sets:

Primer/probe set	Optimal temperature (defined by steepest curve)
# Sarbeco	62.3 - 63
# 1	60.8 - 63
# 5	55.1 - 63
# 6	55.1 - 58
# 7	curves looked good at all temperatures in gradient
# 19	52.6 - 62.3
# 20	52.6 - 58.2

From this, the following temperature gradient was created and included in the thermal (PCR) protocol:


See next page.






# Covid-19 project - Digital lab journal

Temperature gradient:




## Results:

A



B



**Finding the optimal temperature for all primer/probe sets. (A)** Box plot from the experiment. **(B)** Graph showing the regression-fitted curves of all samples (standard deviation included).



### Conclusion:

---

- We will keep testing the primer/probe sets as we did in this project by including the temperature gradient in the thermal (PCR) protocol.
- Optimal temperature did not necessarily match the melting temperature for all primer/probe sets (tested by manufacturer with NaCl).



## Reducing extension time to 30 seconds

SN0010

21.09.2020

### Aim:

Try to reduce extension from 1 minute to 30 seconds (s), to save PCR time.

### Mix:

#### One-step RT-PCR

Fragment:	Control #1	Date	21.09.2020
PCR volum, µL	20	Temp. Gradient 55-62°C	# of reactions 13


Reagents	Working sol. conc.		Total volume 260 µL		Desired conc.	
			Volume			
H <sub>2</sub> O			208.52			
10X Thermopol uten MgSO <sub>4</sub>	0		26.00			
MgSO <sub>4</sub>	200	mM	2.60		2	mM
Primer forward	100	µM	1.04	Or 3.12 µL from mix	0.4	µM
Primer reverse	100	µM	1.04		0.4	µM
SYBR green	100	µM	1.04		0.4	µM
dNTP	100	mM	1.04		400	µM
RNA/cDNA	10	ng	13.00	12.00	0.5	ng/µl
BSA	100	%	2.60		1	%
Dave Taq native	37.5	U/µl	0.52		0.075	U/µl
Mashup	7.5	U/µl	2.60		0.075	U/µl

### PCR protocol:

- Same as in SN0009


### Results:

A



## Covid-19 project - Digital lab journal

**B**



**Reducing extension time to 30 seconds.** (A) The graph shows the new set-up for all primer/probe sets, including the reduced extension time. (B) Box plot of the experiment.

### Conclusion:

- The set-up optimized so far will be called "Set-up 10" and includes:
  - Homemade kit
  - Temperature gradient when using different primer/probe sets
  - One-step RT-PCR thermal protocol with 5 min cDNA synthesis time and 30 seconds extension time

