

Test of 3 different Taq enzyme in Sars\_CoV 2 assay.

Starting concentration in the first well, serially diluted with the same mixture without TAQ. Thus, the enzyme concentration I halved in each following reaction.

Well # Microliter enzyme/20 microliter reaction mix


1	1,0000
2	0,5000
3	0,2500
4	0,1250
5	0,0625
6	0,0313
7	0,0156
8	0,0078
9	0,0039
10	0,0020
11	0,0010
12	0,0005

Fragment: Sars CoV 2	PCR			Date	16.06.2020	
PCR volum, µl	20		# of reactions		1	
	Working					
	solutions		Total volume 20 µl		Desired	
	concentration		Volum		concentration	
H2O			14,48			
10X Thermopol	20 mM	MgCl	2,00			
MgSO4	200 mM		0,00		2,00001 mM	
Primer Forward	100 µM		0,08		0,4 µM	
Primer reverse	100 µM		0,08		0,4 µM	
Probe	100 µM		0,08		0,4 µM	
dNTP	100 mM		0,08		400 µM	
cDNA	10 ng		2,00		1 ng/µl	
BSA	100 %		0,20		1 %	
Dave/Swiss/Hot -TAQ	20 U/µl		1,00		1 U/µl	
Mashup	7,5 U/µl		0,00		0 U/µl	

Primer and probe

"RdRp gene/nCov\_IP2"  
nCoV\_IP2-12669Fw ATGAGCTTAGTCCTGTTG  
nCoV\_IP2-12759Rv CTCCTTTGTTGTGTTGT  
nCoV\_IP2-12696bProbe(+) /56-FAM/AGATGTCTTGTGCTGCCGGTA/3BHQ\_1/

Dave Taq



Swiss Taq

